ProTime

Microcoagulation System

PROTIME CUVETTE AND PROTIME3 CUVETTE

PACKAGE INSERT

ENGLISH

IVD

The ProTime™ Microcoagulation System consists of a portable, battery operated instrument and disposable cuvette for quantitative determination of prothrombin time (PT) from fingerstick whole blood or anticoagulant-free venous whole blood. There are two different types of cuvettes available with the ProTime Microcoagulation System a ProTime cuvette and a ProTime3 cuvette. The ProTime cuvette is black color coded and utilizes three channels for the PT assays and two channels for integral controls. It requires approximately 65 µl of blood (approximately 3 drops). The ProTime3 cuvette is blue color coded and utilizes one channel for PT assay and two channels for integral controls. The ProTime3 cuvette requires approximately 27 µl of blood (approximately 1 large drop). The ProTime Microcoagulation System is designed for professional use in the management of patients treated with oral anticoagulants or for patient self-testing.

SUMMARY AND EXPLANATION

Traditional coagulation tests measure the time required for the formation of a fibrin clot following the addition of a coagulation activating reagent. Laboratory assays typically use plasma recovered from anticoagulated (clotted) whole blood samples. More recently, microcoagulation systems which utilize small samples of whole blood, such as the ProTime Microcoagulation System, have been used clinically. The clotting time is a measure of the functionality of the patient's hemostatic system. Specific coagulation activating reagents are employed in various clotting tests to measure specific components of the hemostatic system. Clotting times are prolonged in the case of either decreased procoagulant activity or increased anticoagulant activity.

The events leading to the formation of a fibrin clot are simplified in coagulation theory into two coagulation pathways: the intrinsic and the extrinsic pathways. There are two tests clotting factors and proteins involved in this cascade scheme, numbered I through XIII, excluding VI. The PT test measures the extrinsic coagulation pathway and is sensitive to coagulation factors VII, X, V, II and Fibrinogen (I). With the exception of Factor V, vitamin K is a required co-factor for biosynthesis of these factors in the liver. The PT test uses thromboplastin as the active reagent to initiate the extrinsic pathway. The PT test will be prolonged in patients with liver disease or vitamin K deficiency. The test is widely used to monitor oral anticoagulant therapy which suppresses the synthesis of vitamin K-dependent clotting factors.

PRINCIPLE OF OPERATION

The ProTime Microcoagulation System measures the PT using fibrin clot formation and detection. ProTime draws a precise volume of blood into each microchamber of the cuvette which contains thromboplastin, stabilizers and buffers. A sensitive thromboplastin was chosen for low ISI (approximately 1.0), which optimizes the sensitivity of the system. A photo-optic array detects the motion of sample/reagent mixtures as they move through a precision restrictor in each channel. The blood is pumped back and forth until a clot begins to form, obstructing the channels and slowing the flow of blood. ProTime detects the clot when the blood movement decays below a predetermined rate.

Both the ProTime and the ProTime3 cuvettes were designed to perform on-board controls with each test. All cuvettes have two levels of on-board controls in addition to the channels for PT assays. The controls ensure that the test procedure was performed correctly and that the reagent system is functioning properly with every test.

CALIBRATION

The ProTime instrument and cuvettes are pre-calibrated. No additional calibration is required.

REAGENTS

ProTime cuvettes are pre-loaded with dried thromboplastin with added stabilizers and buffers. The thromboplastin has a high sensitivity, measured as an ISI near 1.0. Each cuvette performs the PT assay and, in addition, has one channel for a Level I control and one channel for a Level II control. The controls consist of purified plasma-extracted coagulation factors and anticoagulants.

PRECAUTIONS

- The ProTime instrument is designed for use only with ProTime cuvettes. ProTime will not produce a result if cuvettes are past their expiration date.
- DO NOT expose ProTime to extreme temperatures above 35°C (95°F).
- Patient specimens and used cuvettes are potentially infectious. The cuvettes include materials that have been prepared from human plasma or serum which has been tested using US FDA licensed methods and found to be non-reactive for HIV antibody and for hepatitis B surface antigen. Handle with appropriate care and disposal of cuvettes and blood collection materials in accordance with standard methods of biohazard control.
- As with all diagnostic tests, ProTime Microcoagulation System test results should be scrutinized in light of a specific patient's condition and anticoagulant therapy. Any results exhibiting inconsistency with the patient's clinical status should be repeated or supplemented with additional test data.
- Do not use if pouch is damaged or open.
- Dispose after use. The ProTime cuvette is for single use only and are not to be reused. After use it contains human blood and should be disposed of in accordance with local regulations for human blood contaminated waste.

REAGENT PREPARATION AND STORAGE

- ProTime cuvettes are ready to use. No additional preparation is required.
 - Store the foil-pouched cuvettes refrigerated 2 - 8°C (36 - 46°F).
 - When stored at 2 - 8°C (36 - 46°F), an unopened cuvette is stable until the date printed on the pouch. Unopened cuvettes may be stored at room temperature (15 - 30°C, 59 - 86°F) for 60 days. Once the pouch has been opened, the cuvette must be used within 16 hours.
- SPECIMEN COLLECTION**
- Fingerstick whole blood is the recommended specimen. Fingerstick and blood collection devices are provided. The Tenderlett™ Plus device is supplied with the ProTime cuvette and the Tenderlett Plus LV (low volume) device is supplied with the ProTime3 cuvette. The Tenderlett Plus device will collect approximately 65 µl of blood (approximately 3 drops) while the Tenderlett Plus LV device will hold approximately 27 µl of blood (approximately 1 large drop). Samples should be analyzed immediately after collection. No additional sample preparation is required.

For venous samples, collect venous whole blood into an anticoagulant-free plastic syringe in place of fingerstick sampling steps 3 and 4 below. Immediately dispense venous sample into the Tenderlett Plus collection cup. Fill the Tenderlett Plus collection cup to the line or the Tenderlett LV collection cup to the top. Follow instructions from step 5 of the fingerstick test procedure.

NOTE: Serum, plasma or whole blood collected with any anticoagulant are NOT suitable samples.

MATERIALS PROVIDED

- ProTime cuvettes with Tenderlett Plus, or
- ProTime3 cuvettes with Tenderlett Plus LV
- Product Instructions
- Antiseptic gauze (with self-test package only)

Information for Patient Self-Testers

For a more detailed description of the device, please refer to your ProTime Operator's Manual. Follow your doctor's instructions if you have difficulty performing the test or if you receive a result outside your therapeutic range.

TEST PROCEDURE (with fingerstick sample)

- Turn on ProTime**
 - Depending on your instrument model, press either the **●** or **○** button. ProTime performs a self-check procedure that may take up to 60 seconds. ProTime will prompt you through the test. Watch the screen and follow the prompts.
- Insert a cuvette**
 - Make sure the ProTime cuvette is brought to room temperature before use. Wait for the prompt. Insert the cuvette into the slot with the printed side face up and the bar code down. The Display will show **WARNING**.
- Prepare for finger incision**
 - Wash your hands with soap and water. Prepare the finger. Wait for the prompt before incising the finger and collecting blood.
 - It is easier to collect blood of the hands are warm. Follow these steps to ensure a good sample:
 - Wash hands in warm water, or rub hands together to stimulate blood flow.
 - Apply firm pressure to the palm and finger. Massage the hand and push blood into the fingertips.
 - Cleanse the middle or ring finger and dry. To prevent contamination, do not touch the site after clotting.
- Blood collection**
 - Wait for the prompt before incising the finger.
 - Place the Tenderlett Plus firmly against the side of the finger and press the red trigger.
 - Wipe away the first trace of blood. Gently massage from the base of the finger to force blood to the tip so that a large drop of blood forms.
 - Touch the large drop of blood to the collection cup. Keep adding blood until the blood level fills the cup above the line.
 - For Tenderlett Plus LV, ensure the cup is filled completely. Ensure the blood extends all the way to the bottom of the cup and either drop it if you are not sure you have enough.
- Snap Tenderlett Plus to ProTime**
 - Hold the device at an angle and place the front end of the device into the slot in the instrument.
 - Press down to click the Tenderlett Plus in place. You should hear a soft click.
- Start the test**
 - Depending on your instrument model, press either the **●** or **○** button. This signals ProTime to draw the sample into the cuvette. It takes only a few seconds for ProTime to draw the blood into the cuvette. Watch the screen for the next prompts.
- Remove Tenderlett Plus**
 - Remove Tenderlett Plus immediately when prompted to do so.

CAUTION: Failure to do so will result in an error message. ProTime allows you six seconds.

- Do not press the **●** or **○** button after the Tenderlett Plus is removed from the ProTime while you are testing. This will interrupt the test procedure, and you will have to start over with a new blood sample.
 - The instrument then progresses to the test and displays **TESTING** screen.
- 8. Read the result.**
- Read the result as displayed on the screen.
 - Depending on your instrument model, press either the **●** or **○** to turn off the instrument.
 - Press either the **▼** or **□** button to go to the **MAIN MENU** if you want to run another test, review the data in memory, print results, transfer results to a computer, or perform set up functions.

NOTE: ProTime measures normal and therapeutic PT's in fresh whole blood. Results are displayed in plasma equivalent times and INR.

CONVERSION EQUATION

The programmed equations used to convert the whole blood clotting time to INR were determined from multi-site clinical trials using a split sample design for both the standard ProTime and the ProTime3. In these studies, whole blood samples were tested immediately using the ProTime system. Plasma samples were prepared from citrated whole blood and were tested in a reference laboratory using conventional laboratory instrumentation and reagents. The conversion equations were obtained by regression analysis comparing the whole blood clotting time to the reference lab plasma result.

The standard ProTime conversion equation is correlated to laboratory results using plasma obtained with tubes containing 3.8% sodium citrate. The ProTime3 conversion equation is correlated to laboratory results using plasma obtained from tubes containing 3.2% sodium citrate. Clinical results demonstrate a high correlation between the standard ProTime and the ProTime3. Slight differences in INR values between the ProTime products are in part attributable to the known difference in plasma INR results with 3.2% or 3.8% sodium citrate.

NOTE: For further information regarding the effect of sodium citrate refer to Suggested Reading, article by Adcock et al.

INTERNATIONAL NORMALIZED RATIO

In order to standardize for the variable potency of different thromboplastins, the International Council for Standardization in Hematology and the International Committee for Thrombosis and Hemostasis have proposed that each manufacturer of thromboplastin reagent identify the potency of their preparation compared to the International World Health Organization (WHO) reference (ISI = 1.0). Such comparison is the basis for assigning the International Sensitivity Index (ISI) of the reagent. International Normalized Ratios (INRs) are calculated using the equation INR = (PT/PT_{ref})^{ISI}, where PT = prothrombin time, PT_{ref} = normal population PT at the local lab and ISI = potency of the plasma reagent system. The ProTime INR is calculated from the whole blood clotting time and the plasma equivalent PT = calculated from the INR (using ISI = 1.0).

The ProTime Microcoagulation System uses a thromboplastin with a low ISI (approximately 1.0) which optimizes the sensitivity of the system.

SENSITIVITY

ProTime is sensitive to deficiencies in vitamin K-dependent coagulation factors known to influence the PT test (ie., Factors I, V and X).

Hematocrit levels between 20% to 60% do not significantly affect test results.

REPORTABLE RANGE

In clinical trials, no significant difference was observed between fingerstick and venous specimens run on ProTime. ProTime measured patients with an INR range of 0.8 to 7.0. ProTime is programmed to calculate INR from 0.8 to 9.0. If INR > 7.0, the numerical result is marked with an "x". If INR > 9.0 or < 0.8, a message indicating this is displayed.

PERFORMANCE CHARACTERISTICS

Precision

Precision testing was conducted with two levels of standard control plasma substrate preparations (Level I and Level II).

Control	Precision	n	Mean	SD
Level I	Within Day	17	0.9	0.06
	Day-to-Day (5 days)	4/day	1.0	0.08
Level II	Within Day	19	3.2	0.19
	Day-to-Day (5 days)	4/day	3.2	0.12

ProTime3 Cuvette	Precision	n	Mean	SD
Level I	Within Day	18	0.9	0.07
	Day-to-Day (5 days)	4/day	0.9	0.12
Level II	Within Day	20	4.0	0.19
	Day-to-Day (5 days)	4/day	4.2	0.22

Accuracy

INR results generated by the ProTime and ProTime3 cuvettes using venous and fingerstick whole blood samples were compared to INR values obtained using standard Laboratory Plasma PT methods with samples collected in 3.2% sodium citrate tubes. The following accuracy data were obtained.

Blood Sample	Linear Equation	r	n
fingerstick	y = 0.94x + 0.38	0.95	229
venous	y = 0.91x + 0.44	0.94	232

ProTime3 Cuvette vs. Lab (Plasma)	Linear Equation	r	n
Blood Sample			
fingerstick	y = 1.05x + 0.07	0.95	229
venous	y = 0.97x + 0.19	0.95	219

Clinical Performance Comparison
In a trial of patient self-testing (PST) in the home vs. professional testing in the clinic and reference lab, equivalent ProTime results were obtained.

ProTime professional fingerstick vs. patient self-testing

y = 0.94x + 0.13 r = 0.92 n = 421

ProTime patient self-testing vs. reference lab

y = 0.77x + 0.38 r = 0.87 n = 368

ProTime3 vs. ProTime
Linear regression containing clinical fingerstick results from three sites yielded a regression equation as below:

y = 1.02x - 0.14 r = 0.94 n = 229

NOTE: For further performance characteristics refer to Suggested Readings, studies published by the Oral Anticoagulation Monitoring Study Group.

QUALITY CONTROL

The ProTime instrument is designed with multiple systems to ensure proper instrument function. The instrument self-check is actuated by check temperature and timing functions, battery level, and optical, electrical and mechanical functions. The instrument does not require further calibration. Each ProTime cuvette has two integral reagent controls which ensure assay reliability and performance. Both levels of control produce quantifiable clotting endpoints which are compared to pre-set acceptance limits programmed on the cuvette barcode. Other in-process instrument QC features and the integral reagent controls function together to ensure correct sample size, collection technique, test procedure, instrument functionality and reagent integrity. When an error message is displayed, the user should review the product instructions and repeat the test.

Additional external control materials may be used to check instrument function, reagent integrity and user technique. Each institution should establish a quality assurance program for PT testing which complies with existing local, state and federal regulations as applicable.

LIMITATIONS

- ProTime uses only fresh capillary or venous whole blood. Plasma or serum cannot be used.
- Glass tubes or syringes must not be used. Only only plastic syringes without anticoagulants to collect venous samples.
- Reagents may be limited in patients receiving supra-therapeutic heparin or who have an abnormal response to heparin.
- Poor fingerstick blood collection technique may affect results.
- Correlation of results reported by the ProTime Microcoagulation System to local laboratory results depends on the precision of the laboratory method and the specific laboratory reagent and instrument system employed.
- Studies indicating interferences may occur in patients with antiphospholipid antibodies or antiphospholipid syndrome.
 - Sanfilippo MJ, Sennet J, McMahon EJ. 2000. Falsely elevated INRs in warfarin-treated patients with the lupus anticoagulant. *WMJ*. Jun; 99(3):62-4.
 - Rosborough TK, Shepherd MF. 2004. Unreliability of international normalized ratio for monitoring warfarin therapy in patients with lupus anticoagulant. *Pharmacotherapy*. Jul; 24(7):838-42.
 - Perry SL, Samsa GP, Orth TL. 2005. Rate of false case testing of the international normalized ratio in patients with antiphospholipid antibodies. *Thromb Haemost.* Dec; 94(6):1196-202.

More Information
For further information and technical assistance, contact ITC Technical Support at 1-732-548-5700, 1-800-631-5945 or techsupport@tmed.com.

Sanfilippo MJ, Sennet J, McMahon EJ. 2000. Falsely elevated INRs in warfarin-treated patients with the lupus anticoagulant. *WMJ*. Jun; 99(3):62-4.

Rosborough TK, Shepherd MF. 2004. Unreliability of international normalized ratio for monitoring warfarin therapy in patients with lupus anticoagulant. *Pharmacotherapy*. Jul; 24(7):838-42.

Perry SL, Samsa GP, Orth TL. 2005. Rate of false case testing of the international normalized ratio in patients with antiphospholipid antibodies. *Thromb Haemost.* Dec; 94(6):1196-202.

EN SUGGESTED READING
DE LITERATURHINWEISE
ES LECTURAS SUGERIDAS
FR LECTURES CONSEILLÉES
IT BIBLIOGRAFIA
NL AANBEVELING LEESMATERIAAL

- Adcock, D.M., Keesin, D. and Marlar, R.A. *Clin Pathol* 1997; 107:103-110.
- Brien, WF., Baskerville, J.C., Taberner, D.A. and Crawford, L. *Am J Clin Pathol* 1999; 111: 199-201.
- Eckman, M.H., Levine, H.J., Paolker, S.G. *Chest* 1992; 102:585-5485.
- Eckman, M.H., H. Levine, S.G., Paolker, N.E. *NEJ* 1993; 329:96-102.
- Fairweather, R.B., Ansell, V., van den Beselae, A.M.H.P., et al. *Arch Pathol Lab Med* 1998; 122:678-781.
- Gosselin, R., Owings, J.T., White, R.H., et al. *Thromb Haemost* 2000; 83:698-703.
- Hirsch, J., Dalen, J.E., Deykin, D., Poller, L. *Chest* 1992; 102:1235-3265.
- Hirsch, J., Poller, L. *Arch Intern Med* 1994; 154:282-288.
- Hubbard, R.R., et al. *Brit J Haemat* 1999; 104:655-660.
- Hirsch, M.K., Hirsch, J., Landefeld, S., Rasbati, G. *Chest* 1992; 102:3255-3035.
- Oral Anticoagulation Monitoring Study Group. *Am J Clin Pathol* 2001; 115:280-287.
- Oral Anticoagulation Monitoring Study Group. *Am J Clin Pathol* 2001; 115:288-296.
- Tripathi, A., Chantarangkul, V., Bressi, C., Manucci, P. *Am J Clin Pathol* 2003; 119:241-245.
- White, R.H., Zhou, H., Ramano, P., Mungall, D. *Clin Pharmacol Ther* 1995; 58:588-93.

8 Olsen Avenue • Edison, NJ 08820 USA
Tel: 732 548 5700 • fax: 732 248 1928

www.tmed.com
a subsidiary of Thoratec Corporation

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